



CONTENT

I. USER MANUAL

1. Technical information	3
1.1 Package content	3
1.2 Technical Specifications	3
1.3 Main components	4
1.4 Installation	4
1.5 Operation rules	6
1.6 Sterilization/Disinfection	7
1.7 Regular maintenance	7

II. WARRANTY

2. Warranty policy	11
2.1 Warranty scope	11
2.2 Manufacturer's warranty	12
2.3 Terms of delivery	13
2.4 Disposal information	14
2.5 Administration offices	14

1. Technical information

1.1 Package content

1.	The main unit PR-2000	1 pc.
2.	Power cors	1 pc.
3.	Pressotherapy suit (48 channels)	1 pc.
4.	Inflatable tubes	5 pcs.
5.	User manual	1 pc.
6.	Transport box	1 pc.
Consumables		
7.	Fuse	2 pcs.

1.2 Technical Specifications

Function of the device	pressotherapy
Pressure power	120-240 mm Hg
Pump time	6-10 sec
Treatment areas	48-channels(48 airbags)
Programs of pressotherapy	5 pre-programed modes
Control panel	push-button interface
<i>Network parameters</i>	
Power supply	220-240 V / 50 Hz; 100-120 V / 60 Hz
Power input	115 W
<i>Safety</i>	
Protection class	IP31
<i>Dimensions</i>	
Net dimensions(W*L*H)	470*610*1050 mm
Net weight	31 kg

1.3 Main Components



The main components of the machine as it is shown above:

1. The main unit PR-2000;
2. Pressotherapy suit (48 air channel);
3. Inflatable tubes;
4. Power cord.

1.4 Installation

To install the PR-2000 machine, proceed as follows:

1. Connect the machine to the mains using the power cord.
2. Connect the inflatable tubes of the pressotherapy suit to the connectors on the suit cuffs with the corresponding numbers. To connect the air supply system from the device to the pressotherapy suit, use the connecting cables of the inflatable tubes. Cables 1; 2; 3; 4; 5 are connected to the connectors 1; 2; 3; 4; 5 on the rear panel of the device, and cables A, B, C, D, E, F, G, H, I, J are connected to connectors A, B, C, D, E, F, G, H, I, J of the connecting cables of the pressotherapy suit.
3. Make sure that all the connections are firmly in place and the suit is airtight (does not allow air).
4. Power up the device by pressing the Power switch on the back of the device
5. Choose one of three areas:
 - Arms;
 - Waist;
 - Legs.



6. Select one of five programs on the device's control panel interface using the "SELECT MODE" button: improve blood circulation - gentle mode (M1), lymphatic drainage massage (M2), eliminate swelling and lose weight (M3), accelerated weight loss program (M4), lymphatic drainage accelerated massage (M5), which corresponds to the physical condition of the patient.
7. Set treatment time. Use «Time» regulator.
8. Set the intensity level (pumping). Use «Intensity» button on the control panel;
9. Now the machine is ready for operation;
10. Press «Start/Pause» button to start the treatment.

1.5 Operation rules

1. Please connect the pressotherapy suit when the machine is switched off!
2. Connect / disconnect each wire of the device by holding the wire at its base in the place of the connection / disconnection being performed.
3. **It is forbidden to wipe the acrylic parts of the device with alcohol-containing liquids or detergents.**
4. Do not allow the unit and its accessories to fall.
5. Do not install near any heat sources. Do not expose the unit to direct sunlight.
6. Do not let liquids, gels, etc. sip into the handpiece or input connectors of the machine.
7. During operation, keep free tension of the handpiece wires.
Avoid bending the wires.
8. Parts of the device in contact with the patient's body (tip of the treatment handpiece) should be disinfected before and after operation.
9. Perform the treatment only with necessary skin care products, produced for this type of treatment.
10. To carry out procedures with this unit, the patient must always have a disposable treatment suit.



11. It is forbidden to repair and maintain the device outside authorized service centers.
12. After transportation at a low temperature, before starting operation, the device must be kept for at least 3 hours at a temperature of from +10 to +35 ° C.
13. It is forbidden to use working tools that are not included in the package content.
14. Do not operate the device with a damaged casing, cable and accessories.
15. The device is designed to 8 hours of continuous operation per day. In case of exceeding the specified continuous operating time, this may lead to a malfunction of the device (overheating of vacuum pump).
Therefore, after each session, it is necessary to turn off the device for at least 20 minutes and give it a temporary break from work.
16. Do not start the procedure without first cleaning the working accessories.
17. Do not start the procedure with the patient if he has a predisposition to one of the contraindications to the use of the device.
18. Only specially trained personnel is allowed to carry out pressotherapy procedures.

1.6 Sterilization/Disinfection

1. After each treatment, clean the inside surfaces of the pressure therapy suit by wiping them with a non-abrasive cloth moistened in medical 70% alcohol. Those parts that have come into contact with the patient should be cleaned after each patient with standard disinfecting methods for devices that come into contact with skin. It is recommended to use medical alcohol at a concentration of at least 70%.
2. Do not use disinfectants containing chlorine components.

1.7 Regular maintenance

1. When cleaning the surface of the device, wipe it with a soft cotton cloth moistened with water or soapy water. **Do not use solutions and liquids containing components of chlorine or alcohol.**
2. Clean the treatment room regularly, in order to avoid dusting and contamination of the device, which may lead to malfunction of the device.
3. Before operation, make sure that there are no mechanical damages to the wires, working accessories, and the case of the device and correct connection of working accessories.
4. If necessary, calibrate the screen by repeatedly pressing the empty space on the screen when the machine is turned on. As a result, a blue calibration mode window will appear on the screen, in which you must alternately click on the crosses in the corners of the screen. After this, the calibration will be completed.

In case that a malfunction does not relate to the ones considered above or that it was not possible to eliminate one in accordance with the indicated recommendations, be sure to contact the service center.

II. Warranty

2.1 Warranty scope

1. This warranty is valid only if the warranty card is correctly filled in by the seller, in particular, if there is a clear seller stamp and the date of sale. The warranty period of the product is calculated from the date of purchase. The conditions of the guarantee operate within the framework of the Law of the EU "On Protection of Consumer Rights", are governed by the legislation of the country and in no way restrict the rights of consumers.

2. The warranty expires in the following cases:

- use of the device in violation of the requirements of the Operation Manual;
- damage caused by intentional or erroneous actions of the consumer;
- presence of mechanical or other damage to the product;
- penetration of liquid, dust, insects and other foreign objects into the product;
- disassembly or any other intrusion into the design of the device.

3. The warranty does not cover defects (faults) of the product due to the following reasons:

- normal wear and tear of parts with a limited service life, as well as consumables (batteries, cuffs, etc.), damage that was caused by improper or careless use, storage, maintenance of equipment, improper use or improper installation;
- the use of low-quality components that have developed their working life;
- The action of force majeure (accident, fire, flood, mains failure, etc.).

In accordance with the Law "On Protection of Consumer Rights" "... a manufacturer or organization that performs the functions of a manufacturer on the basis of an agreement with him is obliged to accept goods of inadequate quality at the consumer, and, if necessary, conduct a quality check of the goods. The consumer has the right to participate in quality control goods

In the event of a dispute about the causes of deficiencies of the goods, the manufacturer or organization that performs the functions of the manufacturer must examine the goods at its own expense. The consumer has the right to challenge the conclusion of such an examination in court.

If, as a result of the examination of the goods, it was found that the defects arose after the goods were transferred to the consumer due to their violation of the established rules for using, storing or transporting the goods, actions of third parties or force majeure, the consumer must reimburse the manufacturer or organization performing the functions of the manufacturer to bear the costs of the examination also associated with its holding costs for storage and transportation of goods.



2.2 Manufacturer's warranty

Manufacturer's warranty is subject to the terms and conditions of operation, transportation and storage.

1. Warranty 18 months applies only to the internal configuration of the device;
2. Warranty 6 months applies only for accessories.

The product verification protocol for compliance with regulatory requirements is available upon request.

In the event of a malfunction, the user is entitled to free repair of equipment during the warranty period, subject to the rules of operation. The warranty repair period depends on the complexity of the breakdown and can last from 5 to 30 working days.

Delivery to the service center and back is carried out at the expense of the buyer in the post-warranty period, as well as if the case of a breakdown is not warranty (the reasons for non-warranty cases are listed above). During the warranty period, delivery to the service center is at the seller's expense.

Warranty repairs are made at a service center authorized by the seller.

The service center has the right to refuse warranty repair for:

- violation of the rules of storage, installation and operation;
- presence of mechanical damage;
- absence or violation of seals;
- damage caused by deliberate user actions;
- causing damage caused by ingestion of foreign objects, liquids, insects;
- disassembling the device or performing repairs in an unauthorized service center;
- modification of the design;
- Damage caused by transportation;
- if the warranty period is exceeded.

2.3 Terms of delivery

Delivery and shipment from the warranty center is at the expense of the buyer.

2.4 Disposal information

In cosmetology devices there are no components that can adversely affect human health and pollute the environment. After the end of life, they can be disposed of without special disposal methods. Components (microcircuits, resistors, capacitors etc.) that are operational after completion of the service life may be used as spare parts for repairs of such and similar devices.



2.5 Administration offices

Страна:	Адрес:	Телефон:
Чехия	Prague 5, Štefánikova 18/25; postcode: 15000	+420 (228) 880691
Украина	г. Киев, Ковальский пер., 19	+38 (068) 568 83 66 0 (800) 33 98 14
Россия	г. Москва, Рязанский пр-т, д. 10, стр. 2	+7 (800) 250 5459
Латвия	Rīga, Madonas iela 21-85	+371 (2) 612 2001